
BTL-6000 SUPER INDUCTIVE SYSTEM

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SOMMAIRE

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
REPETITIVE PERIPHERAL MAGNETIC STIMULATION AS PAIN MANAGEMENT SOLUTION IN MUSCULOSKELETAL AND NEUROLOGICAL DISORDERS – A PILOT STUDY

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
PAIN REDUCTION DEMONSTRATED ON THE VISUAL ANALOGUE SCALE

IMPROVEMENT DESCRIBED BY 87% OF THE PATIENTS

REPETITIVE PERIPHERAL INDUCTIVE STIMULATION IN MUSCULOSKELETAL PAIN MANAGEMENT – A PILOT STUDY

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Rehabilitation Center Kladruby, Czech Republic;

Charles University, The Faculty of Physical Education and Sport, Czech Republic

Technical University of Liberec, Czech Republic

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PERIPHERAL APPLICATION OF REPETITIVE PULSE MAGNETIC STIMULATION ON JOINT CONTRACTURE FOR MOBILITY RESTORATION

Kouloulas E.

Rehabilitation Unit "Physiatriki";

2nd Neurosurgery Dpt, University of Athens, "Attikon" Hospital, Athens, Greece

Published: International Journal of Physiotherapy, 2016, 3(5), 519-524

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
EFFECTS OF PARA-SPINAL REPETITIVE MAGNETIC STIMULATION ON MULTIPLE SCLEROSIS RELATED SPASTICITY

Serag H., Abdelgawad D., Emara T., Moustafa R., El-Nahas N., Haroun M.


Neurology Department, Ain Shams University, Cairo, Egypt

Published: International Journal of Physical Medicine and Rehabilitation, 2014, 2:242

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
LUMBAR REPETITIVE MAGNETIC STIMULATION REDUCES SPASTIC TONE INCREASE OF THE LOWER LIMBS

Krause P., Edrich T., Straube A.

Department of Neurology, University of Munich, Klinikum Grosshadern, Munich, Germany

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
REPETITIVE PERIPHERAL INDUCTIVE STIMULATION IN COMPREHENSIVE THERAPEUTIC APPROACH – A CASE STUDY

Žarković D.

Charles University in Prague, Faculty of Physical Education and Sport, Czech republic

Presented : XXIII. Congress of Society for rehabilitation and physical medicine, May, 2016

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FUNCTIONAL MAGNETIC STIMULATION FOR CONDITIONING OF EXPIRATORY MUSCLES IN PATIENTS WITH SPINAL CORD INJURY

Lin V.W., Hsiao I.N., Zhu E., Perakash I.,

The Functional Magnetic Stimulation Laboratory of the SCI Health Care Group, VA Long Beach Health Care System Department of Veterans Affairs, Palo Alto Health Care System

Published: *Archives of Physical Medicine and Rehabilitation*, 2001, 82(2), 162-166

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
RENFORCEMENT MUSCULAIRE

INTENSITY SETTING

Keep verbal contact with the patient while setting the intensity to get feedback on their perception of the contraction.

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

SOLLICITATION ARTICULAIRE

INTENSITY SETTING


Get the patient's feedback while setting the intensity in order to define the level properly.

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FRACTURE

After setting the intensity, the therapy continues operator-free. The therapy is audible while running.

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4

TITLE: REPETITIVE PERIPHERAL INDUCTIVE STIMULATION IN MUSCULOSKELETAL PAIN MANAGEMENT – A PILOT STUDY

Authors: Pěšíký J¹, Váňa Z², Šubert D¹, Zárkovský D¹, Proszka D¹, Bítman V³
 Affiliations: ¹Rehabilitation Center Kladubý, Kladubý, Czech republic;
²Charles University in Prague, Faculty of Physical Education and Sports, Prague, Czech republic;
³Technical University of Liberec, Liberec, Czech republic

ABSTRACT:

Background: Repetitive peripheral inductive stimulation is a solution in musculoskeletal pain management in medical branches such as neurology, orthopaedics, rehabilitation and physical medicine. Therapy is based on the principle of time-varying electromagnetic field passing through neural and muscular tissue, in which electric currents are induced by alternating conductive tissue, a chain reaction including changes in action potential and leading to muscle contraction is performed.

Aim: Aim of the pilot study was to investigate immediate pain relief effect of the repetitive peripheral inductive stimulation device BITL-6000 Super Inductive System (BITL, Industries Ltd.) in musculoskeletal diseases.

Methods: 11 subjects from Rehabilitation Center Kladubý were comprised in the pilot study. Subjects underwent approx. 7 therapies individually. A hand-held applicator type „focus field“ was used. Pulse repetition rates matching with gate theory and peripheral pattern theory were applied to achieve pain relief effect. A Visual Analog Scale (VAS) was used to evaluate pain before and after each therapy. Results: Immediate pain relief effect of the repetitive peripheral inductive stimulation in 62% patients with musculoskeletal diseases was observed.

Conclusion: Despite small number of subjects, decrease of painful perception and pain relief effect were achieved. Although, the study design does not allow comparison with placebo effect, statistically significant ($\alpha = 0.2$) immediate pain relief effect in 62 – 74 % of subjects might be expected.

Key words: Pain, pain relief effect, gate control theory, peripheral pattern theory, repetitive peripheral inductive stimulation, vertebralogenic algic syndrome, visual analog scale

5

TITLE: CLINICAL STUDY ON APPLICATION OF SUPER INDUCTIVE ELECTROMAGNETIC FIELD ON PAIN CONDITIONS

Authors: Štejný F¹, Proszka D¹
 Affiliations: ¹Ortopedie Štátní úst. Praha 6 – Řepy, Child and adult orthopaedic and traumatology clinic, Faculty hospital Motol, Prague, Czech republic;
²Charles University in Prague, Faculty of Sport and Physical Education, Prague, Czech republic

ABSTRACT:

Background: A new approach to treatment of musculoskeletal disorders has appeared in physical therapy. The technology is based on the effect of strong pulse electromagnetic field in human tissue (the induction is measured in the units of Tesla). This pilot study examines the analgesic effect of this technology with different degrees.

Objective: Verification of the analgesic effect of a strong pulse electromagnetic field on a sufficient statistical sample in a clinical practice.

Methods: The therapy was performed with 57 randomly selected patients with chronic and acute pain of musculoskeletal system. Patients had 6 therapies on average, 3-7 times per week, 10-15 minutes according to the selected protocol. We used the combination of the Visual Analog Scale (VAS) and the Verbal Numerical Rating Scale (VNRS) to determine the analgesic effect. The super inductive stimulator (the BITL-6000 Super Inductive System, by BITL, Industries Ltd.) was used to perform the therapy.

Results: Regardless of diagnosis, the overall decrease of pain was 37.5 %. There was significant pain relief by 46 patients. There was neither improvement nor worsening of pain by 4 out of the total 50 patients. 7 patients were excluded from the study.

Conclusion: We have demonstrated the analgesic effect of a strong pulse electromagnetic field on musculoskeletal pain.

3. MUSCLE STIMULATION

TITLE: REPETITIVE PERIPHERAL INDUCTIVE STIMULATION IN COMPREHENSIVE PHYSIOTHERAPEUTIC APPROACH – A CASE STUDY

Authors: Zárkovský D¹
 Affiliations: ¹Charles University in Prague, Faculty of Sport and Physical Education, Prague, Czech republic

ABSTRACT:

Background: Repetitive peripheral inductive stimulation (RPISt) uses high-intensity electromagnetic field, which passes through neural tissue in which electric current causing change in action potential is induced. As the currents are carrying an electric signal to the muscle, muscle contraction is achieved. RPISt is a method covering multiple therapeutic effects in various medical branches. By affecting neuromuscular tissue, RPISt is indicated in treatment of various conditions of musculoskeletal and neural system.

Aim: Aim of this case study was to integrate RPISt in comprehensive physiotherapeutic approach in treatment of patient with posttraumatic respiratory and musculoskeletal disease. Device Super Inductive System (SIS, BITL, Industries Ltd.) was used to achieve pain relief, myorelaxative and breathing enhancement effect.

Methods: A 29-year-old male patient underwent a 4-week comprehensive therapeutic protocol comprised of total 16 therapies. Spirometry and kinology evaluation were used to compare the patient's condition „before and after“ in therapy. A hand-held applicator type „focused field“ was used to treat structural areas. As the device allows setting wide range of therapeutic parameters, various stimulation frequencies were used to achieve different therapeutic effects.

Results: The 4-week protocol with SIS resulted into improvement of spirometric parameters. Statistically significant changes in SVC, FVC and MVV profile were obtained. A positive effect was also observed in the musculoskeletal system, where numerous painful muscle spasms, leading to scoliotic trunk asymmetry were eliminated.

Conclusion: This case study suggest that RPISt can be effectively integrated in comprehensive physiotherapy treatment and covers various therapeutic effects. Although, the study represents only one case, statistically significant changes in respiratory and musculoskeletal system were observed.

Key words: repetitive peripheral inductive stimulation, spirometry, kinology, physiotherapy

4. ANALGESIC EFFECT

TITLE: REPETITIVE PERIPHERAL MAGNETIC STIMULATION AS PAIN MANAGEMENT SOLUTION IN MUSCULOSKELETAL AND NEUROLOGICAL DISORDERS – A PILOT STUDY

Authors: Karasikova K¹
 Affiliations: ¹Prokop Hospital, Štá. Gen. Rebenek, Šofka, Bulgaria

ABSTRACT:

Background and Objectives: Non-invasive therapeutic approaches without negative side-effects are desirable in pain condition treatment where the mobility limiting factor is also there. Repetitive peripheral magnetic stimulation (rPMS) is considered a promising curative method from different perspectives. Because of wide range of therapeutic effects, therapy is mainly indicated in musculoskeletal and neurological disorders. Aim of this study was to investigate pain relief effect and improving of the difficulties in performing Activities of Daily Living (ADL), achieved by rPMS among patients with acute and chronic conditions associated with musculoskeletal and neurological painful disorders.

Methods: 40 patients (n = 23 women, n = 17 men) with acute and chronic painful condition and difficulty to perform ADL, accompanying musculoskeletal or neurological disorders were comprised in the study. All patients were treated with rPMS. The therapy parameters were adjusted to patient's condition. Patients with acute pain underwent daily treatments (n = 5). Patients with chronic pain underwent treatments three times per week (n = 10). The pain presence was evaluated by a 10-point Visual Analog Scale (VAS) for Pain (Pain). Difficulties to perform ADL were evaluated by Patient Functional Assessment Questionnaire (PFAQ) for ADL. A three-month follow-up was completed. All collected data were further evaluated.

Results: There was a statistically significant difference in the before/after condition comparison. Majority of participants described pain decrease (89.33 %) on VAS for Pain and improvement (43.83 %) in ability to perform ADL, after the course of treatment. A three-month follow-up showed bearing improvement (to 43.04 % (vs. before treatment condition)) in ADL performing abilities.

Conclusions: Similar results proved that rPMS therapy can be used as an effective and non-invasive treatment of painful condition with ADL limiting factor accompanying musculoskeletal and neurological disorders. Persisting pain relief effect and ameliorating patient quality of life were observed.

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